




atomization



DART with Syringe



DART without Syringe

Soft conical foam tip to facilitate comfortable application



REACH 8 inch Stylet with Syringe



REACH 8 inch Stylet without Syringe



Flexible stylet that can be adjusted



REACH 4 inch Stylet without Syringe

DART intranasal mucosal atomization device

Deliver medications quickly and painlessly with the DART. Atomized medication delivery results in a therapeutic onset of action similar to intravenous delivery - with no pain or delay associated with starting an IV.

Features:

- Contain soft conical-shaped foam tip to facilitate comfortable application and forms a good seal within the nostril
- Measure and deliver with accuracy
- No IV line needed

item #	description	qty
301-2170-20300	Curaplex [®] DART without syringe	25/ca
301-2170-20150	Curaplex [®] DART with 1mL syringe	25/ca
301-2170-20100	Curaplex [®] DART with 3mL syringe	25/ca

DART Reach laryngo-tracheal mucosal atomization device

Deliver medications quickly and effectively for difficult airways with DART Reach. Features a flexible stylet that gives control over the procedure and the ability to adjust the tip based on the patient's anatomy.

Features:

- Flexible stylet that can be adjusted based on the patient's anatomy
- ISO standard luer lock
- Available in adult and pediatric sizes

item #	description	qty
301-REACH600	Curaplex [®] DART Reach (8 inches) Device with 5mL Syringe	25/ca
301-REACH700	Curaplex [®] DART Reach (8 inches) Device without 5mL Syringe	25/ca
301-REACH720	Curaplex [®] DART Reach (4 inches) Device without 5mL Syringe	25/ca