



Protecta-Gel™

Protecta-Gel™ by Curaplex® was designed by experts in respiratory care to minimize discomfort and skin breakdown in patients prone to pressure injuries. It is hypoallergenic so it is safe for use on thin or sensitive skin. Unlike other gels or barriers on the market, Protecta-Gel™ is specifically designed to fit softly and snugly around the bridge of the nose and cheekbone areas, making it convenient for the clinician and comfortable for the patient. Protecta-Gel™ also features a revolutionary design with patented “ribs”, acting as shock absorbers, providing added comfort for full face mask/nasal mask CPAP and NIV users.



- Latex-free and hypoallergenic - made from a medical grade mineral oil
- Super soft barrier between the patient’s skin and their full face mask or nasal mask
- Unique patented “pressure absorbing ribs” increase comfort and decrease soreness, marks, & wounds to the nasal bridge, upper and lower cheek bone areas
- Can be cut or trimmed to customize sizing for nasal masks
- Tapered ends designed to prevent leaks and improve comfort
- Non-adhesive
- Can be used for up to 30 days with proper care

item #	description	qty
301-T1126	Protecta-Gel™ Skin Barrier	10/bx
301-T1126SM	Protecta-Gel™ Skin Barrier, Small	10/bx