




EtCO₂ sampling cannulas

EtCO₂/O₂ sampling cannulas collect the exhaled CO₂ while simultaneously delivering oxygen. These EtCO₂/O₂ cannulas also feature a reflective style connector, compatible with Microstream™, Capnostream™ and Oridion™ capnography monitors.

- Simultaneously deliver oxygen while obtaining CO₂ sampling during spontaneous breathing
- Soft cannula and tubing material help to ensure patient comfort
- 3-channel tubing prevents kinking
- Available in three tubing lengths, each with a Universal O₂ connector



Microstream™, Capnostream™, and Oridion™ are trademarks of Medtronic, related to Product marked by Medtronic.

standard EtCO₂ cannulas

Standard use dual cannulas available in Nasal and Oral/Nasal configuration.




Nasal Cannula



Oral/Nasal Cannula

extended use EtCO₂ cannulas

Extended use cannulas are fitted with a section of Nafion tubing that allows the cannula to function longer and in higher-humidity conditions.




Nasal Cannula



Oral/Nasal Cannula

sampling line with airway adapter (T-piece)

Sampling line with Airway Adapter (T-Piece) can be used on both intubated and non-intubated patients. Available in both standard use and extended use. Can be used on Adult/Pediatric patients with ETT tube size > 4.5 mm.



Standard Use



Extended Use

Item #	Description	Size		Configuration			Connector/ Length of Use		Length of Sample Line		
		Adult	Pediatric	Nasal Cannula	Oral/ Nasal Cannula	Airway Adapter	Orange (Standard)	Yellow (Extended use)	7'	10'	14'
301-5107F-SE	EtCO ₂ /O ₂ Nasal Cannula	✓		✓			●		✓		
301-5110F-SE	EtCO ₂ /O ₂ Nasal Cannula	✓		✓			●			✓	
301-5114F-SE	EtCO ₂ /O ₂ Nasal Cannula	✓		✓			●				✓
301-5207F-SE	EtCO ₂ /O ₂ Nasal Cannula		✓	✓			●		✓		
301-5210F-SE	EtCO ₂ /O ₂ Nasal Cannula		✓	✓			●			✓	
301-5214F-SE	EtCO ₂ /O ₂ Nasal Cannula		✓	✓			●				✓
301-5707F-SE	EtCO ₂ /O ₂ Oral/Nasal Cannula	✓			✓		●		✓		
301-5714F-SE	EtCO ₂ /O ₂ Oral/Nasal Cannula	✓			✓		●				✓
301-5107F-SE-E	Extended EtCO ₂ /O ₂ Nasal Cannula	✓		✓				●	✓		
301-5114F-SE-E	Extended EtCO ₂ /O ₂ Nasal Cannula	✓		✓				●			✓
301-5207F-SE-E	Extended EtCO ₂ /O ₂ Nasal Cannula		✓	✓				●	✓		
301-5214F-SE-E	Extended EtCO ₂ /O ₂ Nasal Cannula		✓	✓				●			✓
301-5707F-SE-E	Extended EtCO ₂ /O ₂ Oral/Nasal Cannula	✓			✓			●	✓		
301-5714F-SE-E	Extended EtCO ₂ /O ₂ Oral/Nasal Cannula	✓			✓			●			✓
301-5107ET	Sampling Line w/ Airway Adapter Intubated	✓	✓			✓	●		✓		
301-5107ET-E	Sampling Line w/ Airway Adapter Intubated	✓	✓			✓		●	✓		
301-5114ET	14FT Sampling Line w/ Airway Adapter Intubated	✓	✓			✓	●				✓

All item numbers: 25 units per box



tri-anim.com | 800.874.2646