

# POM<sup>®</sup> Procedural Oxygen Mask<sup>®</sup>



The Procedural Oxygen Mask™ (POM<sup>®</sup>) is a versatile multi-port capnography mask that allows most types of scopes, probes, and tubes to be inserted through the mask and enter the nose or mouth while maintaining a high percentage of oxygen delivered to the patient. It maximizes oxygen delivery and provides reliable end tidal capnography monitoring for increased patient safety.

**Capnography Port** allows for continuous EtCO<sub>2</sub> monitoring.



**Penetrable Membranes** allow for scopes, probes, and tubes to access the nose and mouth while maximizing FiO<sub>2</sub>.

Item #	Size		Adapter O <sub>2</sub> Concentration (7ft O <sub>2</sub> Line)		EtCO <sub>2</sub> Configuration (10ft Sample Line)			Qty
	Adult	Pediatric	Medium and High (assembly required)	Medium (preassembled)	Microstream™ Compatible	Male to Male	Male to Female	
301-0318LT-MSA	✓	✗	✓	✗	✓	✗	✗	30/box
301-0318LT	✓	✗	✓	✗	✗	✓	✗	30/box
301-0318LT-FM	✓	✗	✓	✗	✗	✗	✓	30/box
301-0319LT	✓	✗	✓	✗	✗	✗	✗	30/box
301-0318LTEZ-MSA	✓	✗	✗	✓	✓	✗	✗	30/box
301-0318LTEZ-MM	✓	✗	✗	✓	✗	✓	✗	30/box
301-0318LTEZ-MF	✓	✗	✗	✓	✗	✗	✓	30/box
301-0318LTEZ	✓	✗	✗	✓	✗	✗	✗	30/box
301-0318LTMSA-PED	✗	✓	✓	✗	✓	✗	✗	30/box
301-0318LTMM-PED	✗	✓	✓	✗	✗	✓	✗	30/box
301-0318LTFM-PED	✗	✓	✓	✗	✗	✗	✓	30/box
301-0318LT-PED	✗	✓	✓	✗	✗	✗	✗	30/box
301-0318LTEZMSA-PED	✗	✓	✗	✓	✓	✗	✗	30/box
301-0318LTEZMM-PED	✗	✓	✗	✓	✗	✓	✗	30/box
301-0318LTEZMF-PED	✗	✓	✗	✓	✗	✗	✓	30/box
301-0318LTEZ-PED	✗	✓	✗	✓	✗	✗	✗	30/box

Mask Dimensions Adult: W = 3.5 in OD / L = 5.8 in OD    Mask Dimensions Ped: W = 2.9 in OD / L = 3.7 in OD    Material: Poly-vinyl chloride (Latex Free!) DEHP / DINP Free



800.874.2646 • www.tri-anim.com

For more information, please contact your Account Manager or call Customer Service at 800-TRI-ANIM (874-2646).