

With AccuPAP[™],
you can treat
atelectasis effectively
and efficiently.
Just follow these
simple steps.



POSITIVE AIRWAY PRESSURE THERAPY HAS NEVER BEEN EASIER!

Adjustable pressure valves

Eliminates the need for a manometer

Patient end with 22mm OD fitting

Can be used with a mask, trach adapter,
or during nebulization

Removeable mouth piece

Allows for patient's preferred
method of treatment delivery

Gas inlet port

Connects via O₂ tubing
(7' included) to wall air
or O₂ flow meter

Pressure manometer port

For optional pressure monitoring

