



Protecta-Gel™ by Curaplex was designed by experts in respiratory care to prevent discomfort and skin breakdown in patients prone to pressure ulcers. Hypoallergenic in nature, it is safe for use on various types of thin or sensitive skin. Unlike other gels or barriers on the market, Protecta-Gel™ is specifically designed to fit softly and snugly around the bridge of the nose and cheekbone areas, making it convenient for the clinician and comfortable for the patient. Protecta-Gel™ also features a revolutionary design with patented “ribs”, acting as shock absorbers, providing added comfort for full face CPAP and NIV users.

- Hypoallergenic and made from a medical grade, mineral oil based gel
- Super soft barrier between the patient’s skin and their full face mask or nasal mask
- Unique patented “pressure absorbing ribs” increase comfort and decrease soreness, marks, & wounds to the nasal bridge, upper and lower cheek bone areas
- One size fits all — can be cut or trimmed to customize sizing for nasal masks
- Tapered ends designed to prevent leaks and improve comfort
- No adhesives necessary
- Can be used for up to 30 days with proper care
- Latex-free



Ordering Information

Protecta-Gel Skin Barrier

301-T1126

Protecta-Gel Skin Barrier, Small

301-T1126SM

For more information, please contact your Account Manager or call Customer Service at 800-TRI-ANIM (874-2646).



Tri-anim®

800.874.2646 • www.tri-anim.com