

## Model 45 Nebulizer

HCPCS:



- Model 45 hand-held spray nebulizer delivers aerosol particles of uniformly consistent density deep into the lungs The break-resistant nebulizer easily disassembles for sterilization
- Used for treatment of asthmatic patients
- Designed for comfortable administration of fine mists of uniform density to reach deep into lung cavity
- Converts liquid (i e antibiotic, oil extracts, etc ) into aerosol
- Used in respirator and mask test fit kits, and bulk shipped as item 4510
- Break-resistant nebulizer that easily disassembles for sterilization by using one of the following methods: autoclaving, boiling or cold sterilization
- Polycarbonate plastic nebulizer may be used with or without mouthpiece (No 25-31-605)
- Complete with mouthpiece, extra jet, and bulb
- Used by research laboratories

### PRODUCTS

### Description

PRODUCTS	Description
45	Hand-Held Spray Nebulizer

### SPECIFICATIONS

### Description

SPECIFICATIONS	Description
MMAD	5 Microns